



Simplicity	
<b>A</b> Minimum swivel radius – with short arm	680mm 644mm
<b>B</b> Minimum wall to seat back	75mm
<b>C</b> Seat depth	370mm
<b>D</b> Minimum folded width to footplate	285mm
<b>E</b> Minimum open width to edge of footplate	580mm
<b>F</b> Armrest width – external	625mm
<b>G</b> Armrest width – internal	480mm
<b>H</b> Seat back height	400mm
<b>I</b> Footplate to seat height – standard seat option – extended seat post	390-490mm 460-560mm
<b>J</b> Minimum footplate height	50mm
Bottom parking space required	660mm <sup>†</sup>
Minimum track intrusion into staircase	160mm
Minimum track intrusion into top landing (typical at 45°)	150mm*
Maximum weight limit – hinge/slidetrack	140kg/22 stone 120kg/19 stone
Speed	0.12m/sec (approx)
Power supply	24v DC battery
Charger	240v
Staircase angle	27-50°

<sup>†</sup>Increases by up to 140mm with zero intrusion kit.  
\*A 0mm track intrusion kit is available.



BUTE

Handicare, 82 First Avenue, Pensnett Estate  
Kingswinford, West Midlands DY6 7FJ, United Kingdom  
Textphone: (180015) 020 7490 8913 ext 05  
Voicephone: (18002) 020 7490 8913 ext 05  
Global Dealer Support: +44 (0) 844 225 3121  
enquiries@handicare.co.uk www.handicare.co.uk

User support: 08700 118282



Please be aware that actual colours may vary slightly from those shown. All information is provided in good faith but does not constitute a guarantee. Minivator Group reserves the right to alter designs and specifications without prior notice.

## Minivator Stairlifts



# The benefits of a Minivator Simplicity Stairlift from Handicare

If you are thinking about installing a stairlift the Simplicity offers you a safe and reliable solution.

Handicare's Minivator Stairlifts have helped thousands of people overcome the problems posed by going up and down their stairs at home.

The Simplicity straight stairlift offers you all the basic functions that you would expect from a stairlift and, if you have a doorway close to the base of the stairs, you can choose from a manually operated hinged track or the unique SlideTrack

that is only available from Handicare.

Supplied in a neutral shade, the wipe clean upholstery will compliment most colour schemes.

The track is fitted to the stairs, not the wall, so it is unlikely you will need to redecorate following the installation of your Simplicity stairlift.

Like all Handicare stairlifts the Simplicity has a battery backup which means it will continue to work in the event of a power cut so you do not need to worry about being stranded up or downstairs.

The compact unit folds up leaving the staircase free for other users and the Simplicity also comes with two remote controls which means, once folded, you can move it along the track freeing up space in your hallway.



## Safety As Standard

Your Simplicity stairlift will come with a range of safety sensors that are designed to stop your stairlift if it meets any obstruction as you travel up the stairs.

Your Simplicity stairlift has a clip operated seat belt for added security and peace of mind.

## Easy to Use

The operating toggle is easy to use, simply hold it towards the direction you wish to travel, up or down, and your Simplicity will move smoothly and gently on its way.

The control is conveniently located at the end of the arm leaving both arms clear to provide unhindered support when getting on or off the seat.

